




Candoo How to iPhone Hearing Impairment & Accessibility

The iPhone makes it easy to stay in the conversation if you're deaf or hard of hearing. There are several accessibility settings which will make your life easier. Made for iPhone hearing aids make it easy to personalize sound settings to your location. Live Listen also helps you have better conversations in loud places.

On the next pages, we will cover a few that we want to highlight.



Hearing Aids Made for iPhone

Apple has worked with top manufacturers to create hearing aids and sound processors designed specifically for iPhone and iPad. These advanced hearing devices provide outstanding sound quality, offer many helpful features, and are as easy to set up and use as any other Bluetooth device.

Here is how you can quickly access the features and settings of your paired hearing device:

- Triple-click the Home button (or the side button on iPhone X or later). Glance at the battery status or change the left and right volume, together or separately.
 - This will quickly apply your audiologist's environmental presets as you go outdoors or enter noisy locations, like restaurants, without having to rely on additional remotes.
 - Click or tap [Here](#) to explore all hearing aid models made for the iPhone.



Noise App (Apple Watch)

Track your hearing health.

The Noise app tracks decibel levels of the ambient sounds around you, helping you identify when the levels could negatively impact your hearing. The information is stored in the Health app so you can refer to your data whenever you need it.







Visible and Vibrating Alerts

See your alerts- Or feel them.


iPhone delivers both visual and vibrating alerts for incoming phone and FaceTime calls, new text messages, new and sent mail, and calendar events. You can set an LED light flash for incoming calls and alerts. You can also have incoming calls display a photo of the caller.




Step 1: Go to Settings



Step 2: Click on accessibility



Step 3: Scroll down to "Audio/Visual"



Step 4: Under LED Flash for Alerts, slide the button on the right to make it green.


Pairing hearing aids with iPhone

Hearing aid support is an accessibility feature that allows the iPhone and iPad to connect with and manage compatible hearing aids.


You can connect to most Bluetooth enabled hearing aids as well as special MFi (made-for-iPhone and iPad) hearing aids.




Step 1: Go into Settings and select Accessibility




Step 2: Select Hearing Devices



Step 3: Make sure your Bluetooth is On



Step 4: iPhone will now search for hearing aids



Step 5: Select the hearing aids when they become available